

S

SURGICAL &
SPECIALTIES

Obstetrics

m|devices®



We believe in effortless performance, so you can focus on caregiving.

One size never fits all. For over a decade we've been designing and manufacturing devices that simplify the roles of clinicians and carers – not just cut to a standard specification but tailored to the unique needs of Australians.

Our focus is on making performance effortless, leaving you free to focus on what matters most – caring for patients.



Leading by Design is a philosophy that reflects our commitment to research and collaboration. A crucial part of our success has been our effort to embed medical professionals in the design process. Their insight helps inform solutions that make the giving of care simple and rewarding.

The problems you face today are the solutions we're working towards for tomorrow.



S

SURGICAL &
SPECIALTIES

Obstetrics


In a clinical context, you need the reliability and performance that comes from solutions carefully tailored to your needs .

Our range of obstetric devices are designed to make the process of giving care simple. Their ease of use enable you to provide a reassuring patient-centred experience.

Cervical Dilating Balloon and Securement Device

Our Cervical Dilating Balloon provides a non-prescriptive alternative for midwives and obstetricians to ripen the cervix, pre-labour.

It's design has been created in collaboration with clinicians and features a single balloon to improve patient comfort and the absence of eyelets creates a closed system to reduce infection risk.



Optional Securement Device
(sold seperately)

| Cervical Dilating Balloon and Securement Device | | | STERILE |
|---|------------------------------------|--------|------------|
| RE-ORDER | DESCRIPTION | VOLUME | CARTON QTY |
| OB010020 | Cervical Dilating Balloon | 80mL | 10 |
| UR060001 | Securement Device (8.5cm x 12.5cm) | | 50 |


- Specifically designed to ripen the cervix
- Single Balloon which can hold up to 80mL
- Avoids pharmaceutically activated hyper-stimulation
- Non-prescriptive alternative which is easy to insert
- Readily removable throughout the process
- Optional Securement Device available
- High contrast position markings
- Batch tracking sticker on packaging
- Latex free

Cervical Dilating Balloon Introducer

The Cervical Dilating Balloon Introducer has been designed in collaboration with clinicians to assist with the insertion of the m|devices Cervical Dilating Balloon.

Its malleable nylon material enables a simple and efficient insertion and the J-loop handle allows the placement of the Cervical Dilating Balloon with minimal movement.

The Introducer can be removed by simply reverse-bending the J-loop handle.



Used with the Cervical Dilating Balloon
(sold seperately)

- Bespoke design to aid the insertion of the Cervical Dilating Balloon
- Malleable nylon material to support insertion
- J-loop handle enables the introducer to be secured in place during use
- Easily manipulated into any required angle or shape
- Designed for easy insertion and removal

| Cervical Dilating Balloon Introducer | | | STERILE |
|--------------------------------------|--------------------------------------|--------|------------|
| RE-ORDER | DESCRIPTION | LENGTH | CARTON QTY |
| OB020000 | Cervical Dilating Balloon Introducer | 43.5cm | 20 |

Cervical Dilating Procedure Pack

Efficiency

Our ‘all-in-one’ pack improves performance by eliminating collection times.

Environmental Impact

Significantly reduces individual packaging waste.

Fit for Purpose

All the necessary items are in one pack, specifically designed to support the insertion of m|devices Cervical Dilating balloon.


- ⌕ Enables faster more efficient setup
- ⌕ Eliminates time spent collecting individual items
- ⌕ Assists with maintaining infection control principles
- ⌕ Significantly reduces individual packaging waste
- ⌕ ‘Fit for purpose’ procedure pack
- ⌕ Inclusion of the m|devices Cervical Dilating Balloon 100% Silicone, closed distal end with high contrast position markings printed on the shaft
- ⌕ Includes batch tracking sticker
- ⌕ Latex free



| Cervical Dilating Balloon Procedure Pack | | STERILE |
|--|--|------------|
| RE-ORDER | DESCRIPTION | CARTON QTY |
| OB010002 | Cervical Dilating Balloon Procedure Pack | 30 |

Cervical Dilating Procedure Pack

| Cervical Dilating Balloon Procedure Pack | | | STERILE | | |
|--|-----------------------------|-----|----------|-------------------------------------|-----|
| QUANTITY | DESCRIPTION | REF | QUANTITY | DESCRIPTION | REF |
| 1 | Sterile Field 120cx120cm | A | 1 | m devices Cervical Dilating Balloon | G |
| 1 | Blue Rectangle Tray | B | 1 | XRD Gauze Swab 10x10cm 8Ply (5) | H |
| 2 | Lubricating Gel Sterile 27g | C | 1 | m devices Securement Device | I |
| 4 | Luer Lock Syringe 20mL | D | 1 | Plastic Forceps Clear | J |
| 2 | Gallipots 120mL | E | 1 | Tracking Label | K |
| 1 | SS Sponge Forceps 25cm | F | | | |



With you every step of the way.

As an Australian, family owned business, we're committed to supporting the best healthcare outcomes for all Australians, with Continence Care, Venous Access and Airway Management solutions.

More than just a manufacturer of devices, we're focused on the end-to-end clinician/patient experience.

From order process, to easy disposal, we do all we can to make our role in the caregiving process seamless.



1800 023 321



m|devices

27 Llewellyn Avenue
Villawood NSW 2163
Australia

T +61 (0)2 8718 2800

F +61 (0)2 8718 2801

mdevices.com.au