

## Disposable Laryngoscopes - Yeescop<sup>TM</sup>

The Yeescop<sup>TM</sup> one-piece, single use laryngoscope incorporates handle and blade eliminating handle and blade incompatibility and avoids the Laryngoscope becoming a source of cross contamination.

- One-piece design avoids handle and blade incompatibility
- Powerful illumination
- Single use design avoids cross contamination
- Simple function testing



| Disposable Laryngoscopes - Yeescop <sup>TM</sup> |                                                  |      |             |            | NON-STERILE |
|--------------------------------------------------|--------------------------------------------------|------|-------------|------------|-------------|
| RE-ORDER                                         | DESCRIPTION                                      | SIZE | HANDLE SIZE | BLADE SIZE | CARTON QTY  |
| AN010000                                         | YEESCOPE <sup>TM</sup> - Disposable Laryngoscope | MAC3 | 138mm       | 139mm      | 10          |
| AN010001                                         | YEESCOPE <sup>TM</sup> - Disposable Laryngoscope | MAC4 | 140mm       | 151mm      | 10          |