

## Insulin Syringes

Our disposable 1 mL and 0.5 mL syringes have fixed needles and are available in two needle lengths: 8mm and 13mm and are offered in 27G, 29G and 30G.

The retractable safety syringe version retracts the needle into the barrel after injection to prevent re-use and needle stick injury.



- ⊗ Orange cap for easy identification
- ⊗ Pre-attached needle
- ⊗ Clear barrel
- ⊗ Bold scale markings

### The retractable syringe:

- ⊗ Permanently disabled needle reduces risk of needle stick injury
- ⊗ Eliminates blood splatter
- ⊗ Audible lock
- ⊗ Snap off plunger
- ⊗ Sharps waste reduction

Insulin Syringes				STERILE
RE-ORDER	DESCRIPTION	VOLUME	IU	CARTON QTY
IV120127	Insulin Syringe with Standard Fixed Needle 27G x 13mm	1mL	100	100
IV120129	Insulin Syringe with Standard Fixed Needle 29G x 13mm	1mL	100	100
IV120130	Insulin Syringe with Standard Fixed Needle 30G x 13mm	1mL	100	100
IV120529	Insulin Syringe with Standard Fixed Needle 29G x 13mm	0.5mL	100	100
IV121127	Insulin Syringe with Standard Fixed Needle 27G x 8mm	1mL	100	100
IV121129	Insulin Syringe with Standard Fixed Needle 29G x 8mm	1mL	100	100
IV121130	Insulin Syringe with Standard Fixed Needle 30G x 8mm	1mL	100	100

Retractable Insulin Syringes				STERILE
RE-ORDER	DESCRIPTION	VOLUME	IU	CARTON QTY
IV158029	Insulin Syringe with Retractable Safety Needle 29G x 8mm	0.5mL	50	100
IV158129	Insulin Syringe with Retractable Safety Needle 29G x 8mm	1mL	100	100
IV150127	Insulin Syringe with Retractable Safety Needle 27G x 13mm	1mL	100	100
IV150129	Insulin Syringe with Retractable Safety Needle 29G x 13mm	1mL	100	100