

Hydrophilic Coated Nelaton Catheters

DESCRIPTION

Our Hydrophilic Coated Nelaton Catheters are made from non-toxic DEHP free and PVC. Their smooth water soluble polymer coating combined with heat polished distal tip supports a smooth insertion without the need of lubricants.

FEATURES & BENEFITS

- Thermo-sensitive PVC
- Heat polished eyelets
- Smooth round tip
- Colour coded connectors
- Available for Paediatric, Female and Male patients


PRODUCT DATA

- Not Made with Natural Rubber Latex
- Not Made with DEHP
- Shelf life: 3 years
- ARTG number: 210992
- Materials: PVC

ISO STANDARDS

- ENISO14155-2:2003 Clinical investigation for medical devices for human subjects—part2:Clinical investigation plans
- ENISO 11737-2:2006 Sterilisation of medical devices-Estimation of the population of micro-organisms on product-Part 2: Guidance
- ISO10993-1:2003 Biological evaluation of the medical devices-Part 1: Guidance on selection of tests

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Tip
2	Side Eyelet
3	Tube
4	Funnel - Colour Coded

Hydrophilic Coated Nelaton Catheters

Hydrophilic Coated Nelaton Catheters						STERILE
RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON/QTY	
UR043001	Hydrophilic Coated Nelaton Catheter Paediatric	30cm	8Fr	Light Blue	25	
UR043002	Hydrophilic Coated Nelaton Catheter Paediatric	30cm	10Fr	Black	25	
UR044000	Hydrophilic Coated Nelaton Catheter Female	18cm	8Fr	Light BLue	25	
UR044001	Hydrophilic Coated Nelaton Catheter Female	18cm	10Fr	Black	25	
UR044002	Hydrophilic Coated Nelaton Catheter Female	18cm	12Fr	White	25	
UR044003	Hydrophilic Coated Nelaton Catheter Female	18cm	14Fr	Green	25	
UR044004	Hydrophilic Coated Nelaton Catheter Female	18cm	16Fr	Orange	25	
UR044005	Hydrophilic Coated Nelaton Catheter Female	18cm	18Fr	Red	25	
UR045000	Hydrophilic Coated Nelaton Catheter Male	40cm	8Fr	Light Blue	25	
UR045001	Hydrophilic Coated Nelaton Catheter Male	40cm	10Fr	Black	25	
UR045002	Hydrophilic Coated Nelaton Catheter Male	40cm	12Fr	White	25	
UR045003	Hydrophilic Coated Nelaton Catheter Male	40cm	14Fr	Green	25	
UR045004	Hydrophilic Coated Nelaton Catheter Male	40cm	16Fr	Orange	25	
UR045005	Hydrophilic Coated Nelaton Catheter Male	40cm	18Fr	Red	25	