

# Hydrophilic Nelaton with Water Sachet Female

## DESCRIPTION

m|devices Hydrophilic Nelaton Catheter with Water Sachet provides an all-in-one product which offers convenience and ease of use, without the need for additional products. The blue holding sleeve offers easy visualisation whilst maintaining aseptic technique. The hydrophilic coating, heat polished closed distal end and bilateral side eyelets support a smooth and painless insertion and efficient drainage of urine.

## FEATURES & BENEFITS

- Hydrophilic coated
- All-in-one product with water sachet
- Blue coloured holding sleeve to maintain sterility during insertion
- Sterile, easy peel packaging with adhesive wall attachment
- Colour coded connectors
- Heat polished eyelets
- Smooth round tip


## PRODUCT DATA

- Female length: 18cm
- Sizes 8Fr-18Fr
- Not Made with Natural Rubber Latex
- Not Made with DEHP
- Material: Nelaton: PVC
- Sachet: Purified water
- Holding sleeve: PE
- Shelf life: 3 years
- ARTG number: 210992

## ISO STANDARDS

- EN ISO 14155-1:2003 Clinical investigation for medical devices for human subjects – part1 General requirements
- EN550:1994 Sterilisation of medical devices – Validation and routine control of ethylene oxide sterilisation
- EN ISO11607-1:2006 Packaging for terminally sterilised medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems

## PRODUCT DIAGRAM



# Hydrophilic Nelaton with Water Sachet Female

Hydrophilic Nelaton with Water Sachet					STERILE
RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON/QTY
<b>UR047000</b>	Hydrophilic Nelaton with Water Sachet - Female	18cm	8Fr	Light Blue	25
<b>UR047001</b>	Hydrophilic Nelaton with Water Sachet - Female	18cm	10Fr	Black	25
<b>UR047002</b>	Hydrophilic Nelaton with Water Sachet - Female	18cm	12Fr	White	25
<b>UR047003</b>	Hydrophilic Nelaton with Water Sachet - Female	18cm	14Fr	Green	25
<b>UR047004</b>	Hydrophilic Nelaton with Water Sachet - Female	18cm	16Fr	Orange	25
<b>UR047005</b>	Hydrophilic Nelaton with Water Sachet - Female	18cm	18Fr	Red	25