

Hydrophilic Nelaton with Water Sachet Male

DESCRIPTION

m|devices Nelaton Catheter with Water Sachet in a male length, also contains a blue coloured holding sleeve for easy visualisation for an aseptic insertion. This all in one product offers convenience and ease of use without the need for additional products. The hydrophilic coating, heat polished closed distal end and bilateral side eyelets support a smooth and painless insertion and efficient drainage of urine.

FEATURES & BENEFITS

- Hydrophilic coated
- All-in-one product with water sachet
- Blue coloured holding sleeve to maintain sterility during insertion
- Colour coded connectors
- Heat polished eyelets
- Smooth round tip


PRODUCT DATA

- Male length: 40cm
- Sizes 8Fr-18Fr
- Not Made with Natural Rubber Latex
- Not Made with DEHP
- Material: Nelaton: PVC
- Sachet: Purified water
- Holding sleeve: PE
- Shelf life: 3 years
- ARTG number: 210992

ISO STANDARDS

- EN ISO 14155-1:2003 Clinical investigation for medical devices for human subjects – part1 General requirements
- EN550:1994 Sterilisation of medical devices – Validation and routine control of ethylene oxide sterilisation
- EN ISO11607-1:2006 Packaging for terminally sterilised medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems

PRODUCT DIAGRAM



Hydrophilic Nelaton with Water Sachet Male

Hydrophilic Nelaton with Water Sachet					STERILE
RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON/QTY
UR048000	Hydrophilic Nelaton with Water Sachet - Male	40cm	8Fr	Light Blue	25
UR048001	Hydrophilic Nelaton with Water Sachet - Male	40cm	10Fr	Black	25
UR048002	Hydrophilic Nelaton with Water Sachet - Male	40cm	12Fr	White	25
UR048003	Hydrophilic Nelaton with Water Sachet - Male	40cm	14Fr	Green	25
UR048004	Hydrophilic Nelaton with Water Sachet - Male	40cm	16Fr	Orange	25
UR048005	Hydrophilic Nelaton with Water Sachet - Male	40cm	18Fr	Red	25