

Intermittent Closed Catheterisation System

DESCRIPTION

Our Intermittent Closed Catheterisation System is sterile and pre-lubricated with Hydrophilic gel, so no additional preparations is needed. The closed system works like a traditional Nelaton Catheter and provides a simple catheterisation that reduces the risk of urinary tract infections.

FEATURES & BENEFITS

- Sterile closed system
- Pre-lubricated
- Heat polished eyelets
- Colour coded
- Smooth round tip


PRODUCT DATA

- Bag size: 1500mL
- Side holes: 2
- Length: 42cm
- Not Made with DEHP
- Not Made with Natural Rubber Latex
- BPA: No
- Material: PVC
- Shelf life: 3 years
- ARTG number: 210992
- Sterile EO

ISO STANDARDS

- ISO 10993-1:2003 Biological evaluation of the medical devices-Part 1: Guidance on selection of tests
- ISO 14155-2:2003 Clinical investigation for medical devices for human subjects—part2:Clinical investigation plans
- ISO 11135-1:2007 Health-care Product Sterilisation Confirmation and Sterilisation by medical ETO verification and routine control

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Cap
2	Tip
3	Interface
4	Gel
5	Nelaton Catheter
6	Bag

Intermittent Closed Catheterisation System

Intermittent Closed Catheterisation System				STERILE
RE-ORDER	DESCRIPTION	GAUGE	BAG CAPACITY	CARTON/QTY
UR050008	Intermittent Closed Catheterisation System	8Fr	1500mL	25
UR050010	Hydrophilic Nelaton with Water Sachet - Coude Tip	10Fr	1500mL	25
UR050012	Hydrophilic Nelaton with Water Sachet - Coude Tip	12Fr	1500mL	25
UR0500014	Hydrophilic Nelaton with Water Sachet - Coude Tip	14Fr	1500mL	25
UR0500016	Hydrophilic Nelaton with Water Sachet - Coude Tip	16Fr	1500mL	25