

Urine Bag 4000mL

DESCRIPTION

m|devices Urine Bag has a large 4000mL volume capacity, non-return valve, 150cm tubing length with bonded stepped connector and integrated hanger with rope.

The urine bag is used to collect and contain urine via a urinary catheter or connected leg bag. This sterile bag is ideal for overnight drainage due to its large volume capacity.

The T-Tap outlet allows for ease of use to empty the drainage bag and the clear incremental measurements allow for visual urine output.

FEATURES & BENEFITS

- 4000mL Large volume capacity
- Long 150cm large Bore Tubing with bonded stepped connector
- Integrated Hanger with Rope for easy attachment to desired location
- Non-return Valve
- Easy to use T-Tap outlet for emptying the bag
- Clear printed incremental markings on bag
- Ideal for overnight drainage


PRODUCT DATA

- Volume: 4000mL
- 150cm tubing length
- Tubing ID/OD: 9 x 1500mm (inlet tubing)
- Material: PVC, PP, PE, Polyester rope, Nylon air filter
- Not Made with Natural Rubber Latex
- DEHP YES = Product contains Phthalates such as the PVC Plasticizer Di-[2-ethylhexyl]-phthalate (DEHP).
- Shelf life: 5 years
- ARTG Number: 267875
- Sterile EO

ISO STANDARDS

- ISO 14971:2012 Medical Device - Applications of Risk Management.
- ISO 10993-1:2009 Biological compatibility test report
- ISO 8669:1996 Urine collection bags

PRODUCT DIAGRAM



ITEM	DESCRIPTION		
1	Cap	6	Air Filter
2	Connector	7	Dog House
3	Tubing	8	Bag
4	Rope	9	Outlet Tubing
5	Hanger	10	T-Tap
		11	Drip Chamber

Urine Bag

STERILE

RE-ORDER	DESCRIPTION	TUBING LENGTH	BAG CAPACITY	CARTON/QTY
UR064000	Urine Bag T-Tap NRV Integrated Hanger with Rope	150cm	4000mL	40pcs